Mölnlycke®

Mepitel® One

WITH SAFETAC® TECHNOLOGY

Single use  
Usage unique

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Soft silicone wound contact layer

## Product description

Mepitel® One is a transparent, perforated and non-absorbent dressing. The perforated structure of Mepitel One allows exudate to pass into an outer absorbent dressing. This can reduce frequent dressing changes closest to the wound and allows for changing the outer dressing with minimised pain and trauma. Mepitel One can also be used as a stand-alone dressing for protection of damaged skin. The dressing has a Safetac® wound contact layer which is a unique adhesive technology. It minimises pain to patients and trauma to wounds and the surrounding skin at dressing removal.

Mepitel One consists of:

- a soft silicone wound contact layer (Safetac) on one side
- a transparent, flexible, thin and perforated polyurethane film

## Intended use

Mepitel One is a wound contact layer designed for the management of a wide range of exuding wounds such as; skin tears, skin abrasions, surgical incisions, partial thickness burns, traumatic wounds, partial and full thickness grafts, radiated skin, leg and foot ulcers.

It can also be used as a protective layer on non-exuding wounds, blisters and on areas with fragile skin.

## Precautions

- If you see signs of infection e.g. fever or the wound or surrounding skin becoming red, warm or swollen, consult a health care professional for appropriate treatment.
- Always consult a health care professional before using Mepitel One on Epidermolysis Bullosa patients.
- When using Mepitel One on partial thickness burns with high risk of rapid granulation or after facial resurfacing, avoid placing pressure on the dressing. Lift and reposition the dressing at least every second day.

- When using Mepitel One for the fixation of skin grafts and protection of blisters, the dressing should not be changed before the fifth day post application.
- Do not use Mepitel One on patients with known sensitivity to silicone or polyurethane.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.

## Instructions for use

Mepitel One can be used by health care professionals or lay-persons.

Mepitel One with this symbol can help you apply and remove the dressing on skin tears. For these wounds, there is a risk of re-opening the flap when removing the dressing. Apply the dressing with the arrow pointing in the same direction as the flap. Start application from the base of the flap.

When removing the dressing, begin removal following the printed symbol, in the direction the arrow is pointing.

1. Cleanse the wound in accordance with clinical practice and dry the surrounding skin thoroughly.
2. Choose a size of Mepitel One that covers the wound and the surrounding skin by at least 2 cm. If needed, the dressing can be cut. For larger wounds, more overlap is required.
3. Remove the protective plastic layers and apply Mepitel One with the sticky side to the wound. The dressing is applied in a correct way when you can read the text printed on the dressing. Smooth the dressing in place onto the surrounding skin to get a good seal. If more than one piece of Mepitel One is used, overlap the dressings. Make sure that the holes are not blocked.
4. Apply an outer absorbent dressing, on top of Mepitel One and fixate.

Minimise overlap and wrinkles, especially when Mepitel One is applied circular or wrapped. This will reduce the risk of stasis of the wound area and limb.

For application guidance see illustrations in this leaflet.

Mepitel One can be used under compression products. Mepitel One can be used in combination with gels.

When the outer absorbent dressing is saturated it should be changed with Mepitel One left in place. Mepitel One is a single use dressing.

Disposal should be handled according to local environmental procedures.

Mepitel One may be left in place for up to 14 days depending on the condition of the wound and surrounding skin, or as indicated by clinical practice. To prevent maceration, exudate should pass freely through the dressing and the holes should not be blocked.

If any serious incident has occurred in relation to Mepitel One, it should be reported to Mölnlycke Health Care.

Mepitel and Safetac are registered trademarks of Mölnlycke Health Care AB.

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Mölnlycke Health Care AG
Mölnlycke Health Care Ltd.
Mölnlycke Health Care US, LLC