


Veress needle

Truly customised laparoscopic solutions




Veress needle is used to create a pneumoperitoneum in laparoscopic procedures. Available both as sterile single-packed product and for inclusion in Mölnlycke® procedure trays.

- Sharp bevelled needle tip ensures simple, safe abdominal wall penetration
- Transparent handle with clear red visual position indicator
- Spring-loaded stylet provides audible and sensory feedback of correct needle placement
- Blunt stylet shields needle tip once correct positioning has been achieved
- High-flow rotating stopcock with luer lock connector for the connection with the insufflation tube
- Available in two lengths: 120mm and 150mm
- Blunt tip with accurately positioned distal insufflation channel for safe and controlled insufflation



Veress needle tip



Veress needle body

| Ref. No. | Description | Size | Pcs/disp. |
|----------|---------------|------|-----------|
| 899105 | Veress needle | 12cm | 25 |
| 899106 | Veress needle | 15cm | 25 |

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00
The Mölnlycke trademark, name and logo is registered globally to one or more of the Mölnlycke Health Care Group of Companies.
© 2017 Mölnlycke Health Care AB. All rights reserved. HQSU0326