

**PRODUCT DATA SHEET**
**MEPIFORM® (NON-STERILE)**

|  |   |   |                        |                                   |
|--|---|---|------------------------|-----------------------------------|
| <b>Product name</b>                            | Mepiform® (non-sterile)   |   |                        |                                   |
| <b>Product description</b>                     | Self-adherent soft silicone dressing for scar care  |   |                        |                                   |
| <b>Intended use</b>                            | Mepiform is designed for the management of both old and new hypertrophic and keloid scars. It can also be used as a prophylactic therapy on closed wounds for prevention of hypertrophic or keloid scarring |   |                        |                                   |
| <b>Product Codes</b>                           | <b>Art.No</b>   | <b>Product Name</b>                               | <b>Size (cm)</b>       | <b>Pack count<br/>CON/RET/TRP</b> |
|  | 293600-00   | Mepiform 4x30 cm                                  | 4x30                   | 1/5/50                            |
|  | 293640-00   | Mepiform 4x30 cm, ph                              | 4x30                   | 1/5/50                            |
|  | 293700-01   | Mepiform 5x7.5 cm                                 | 5x7.5                  | 1/5/25                            |
|  | 293740-01   | Mepiform 5x7.5 cm, ph                             | 5x7.5                  | 1/5/25                            |
|  | 293770-00   | Mepiform 5x7.5 cm, OTC Norden                     | 5x7.5                  | 2/18/144                          |
|  | 293800-01   | Mepiform 10x18 cm                                 | 10x18                  | 1/5/35                            |
|  | 293750-01   | Mepiform 10x18 cm, ph                             | 10x18                  | 1/5/35                            |
| <b>Manufacturer</b>                            | Mölnlycke Health Care AB<br>Box 13080<br>SE-40252 Göteborg<br>Sweden  |   |                        |                                   |
| <b>Country of origin</b>                       | Finland   |   |                        |                                   |
| <b>Product classification according to MDD</b> | Class I   |   |                        |                                   |
| <b>CE-marked</b>                               | Yes, Notified Body BSI 0086   |   |                        |                                   |
| <b>Products material</b>                       | Backing   | Polyurethane film laminated with viscose nonwoven |                        |                                   |
|  | Wound contact layer   | Thin layer of tacky soft Silicone                 |                        |                                   |
|  | Release film  | Polypropylene film                                |                        |                                   |
| <b>Shelf life</b>                              | 3 years   |   |                        |                                   |
| <b>Other info of special interest</b>          | Product and packaging material is free of natural latex   |   |                        |                                   |
| <b>Physical /Chemical characteristics</b>      | Mepiform is released according to the following standards   |   |                        |                                   |
|  | <b>Property</b>   | <b>Reference Method</b>                           | <b>Result</b>          |                                   |
|  | Adhesion to steel   | T-261   | ≥ 0.6 and ≤ 1.6 N/25mm |                                   |

|                          |   |
|--------------------------|---|
| <b>Packaging</b>         | <p>The purpose of the shelf container and transport box is to protect material from dust, moisture and damage during transport, handling and storage.</p> <p>Unit package is a sterile barrier for the product. Do not use if unit package is damaged or opened prior to use.</p> |
| <b>Storage Condition</b> | N/A   |

| <b>GS1 Codes</b> | <b>Art No</b> | <b>CON</b>    | <b>RET</b>    | <b>TRP</b>    |
|------------------|---------------|---------------|---------------|---------------|
|                  | 293600-00     | 7332430998880 | 7332430980199 | 7332430980182 |
|                  | 293640-00     | 7332430998880 | 7332551879754 | 7332551879747 |
|                  | 293700-01     | 7332430998897 | 7323190179411 | 7323190179404 |
|                  | 293740-01     | 7332430998897 | 7323190179206 | 7323190179190 |
|                  | 293750-01     | 7332430998903 | 7323190179237 | 7323190179220 |
|                  | 293770-00     | 7332551949327 | 7332551940959 | 7332551940942 |
|                  | 293800-01     | 7332430998903 | 7323190179268 | 7323190179251 |